

Soldier 2.5



High Flow Microcatheter with Ultra-Thin Wall Technology



Smallest Outer Diameter
(0.027" ID)
2.5 Fr OD

Fastest Flow Rate at 1000 PSI
(2.5 Fr OD)

Exceptional Tracking

High Torquability

Unparalleled Performance


Best in Class Microcatheter

Mastering
Performance
and Delivery

Soldier microcatheters are designed with advanced technology, allowing them to have the highest level of performance of any microcatheter available. Tracking, torquing, curved tips, largest lumen, highest flow rate, high-pressure tolerance...



...get where you want *FAST*, deliver what you want *ALWAYS*.



**Unique 2.5 Fr profile with 0.027" inner diameter
Largest ID, Highest Flow Rate, Best Embolic Compatibility**

Product Code (<u>part number</u>)	Outer Diameter (Distal)	Outer Diameter (Proximal)	Working Length	Inner Diameter	Tip Shape
S27-SP-130	2.5 Fr	2.5 Fr	130 cm	0.027"	Straight
S27-SP-155	2.5 Fr	2.5 Fr	155 cm	0.027"	Straight
S27-KP-130	2.5 Fr	2.5 Fr	130 cm	0.027"	K-Tip
S27-KP-155	2.5 Fr	2.5 Fr	155 cm	0.027"	K-Tip



INDICATIONS FOR USE

The Soldier is indicated for use in the blood vessels of the peripheral vasculature. It is intended to assist in the delivery of diagnostic agents and therapeutic agents into the target treatment area.

CONTRAINDICATIONS

Not intended for use in embolectomy or angioplasty procedures. Not intended for use in the coronary vessels. Not intended for use in the neurovasculature.

Embolx, Inc.
530 Lakeside Dr. #200
Sunnyvale, CA 94085
www.embolx.com

To order product or for more information contact us at customer@embolx.com or +1-408-990-2949.

©2021 Embolx, Inc. All rights reserved.
Soldier is a registered trademark of Embolx, Inc.